




Forensic Fingerprint / Climate Chamber  
for the development of fingerprints  
on porous surfaces, using Ninhydrin,  
DFO and Indaione

Only authorised supplier of Attestor products in Australia  
VISIT [WWW.OPIRA.COM.AU](http://WWW.OPIRA.COM.AU)

## Air Filter and External Ventilation

**Air Filter and External Ventilation** In order to reduce the risk of a contamination of evidence with particulates like dust loaded with chemicals from former treatment cycles, NINcha features a special air filter system. The chamber is also fitted with process controlled flaps for inflow and exhaust for optional connection to external chemical ducts



### Advantages at a glance

- **Pre-Programmed Treatment Cycles** Ready-to-run cycles for Ninhydrin, DFO and Indandione ensure simple operation and consistently high-quality evidence development. Parameters and timers can also be adjusted by a user to suit their needs, or cycles can be controlled entirely manually if desired.
- **Modular Shelving Concept** The brackets on the inner walls of the chamber are designed for the use with both the easy clean stainless steel metal grills and the hanging rods.
- **Air Filter and External Ventilation** In order to limit cross-contamination and simplify cleaning, a special air filter system has been integrated into the chamber to filter out residual chemical developer materials or small evidence particulates that may remain in the chamber after use. The chamber can also be upgraded to connect to an external ventilation system with process-controlled air inflow and exhaust valves.
- **UV-Decontamination Unit** For easy elimination of DNA inside the chamber, a short-wave UV decontamination unit can be installed in the chamber.
- **Water Supply and Drainage** The chamber has an integrated reservoir for distilled water, as well as an automated condensation control system to remove excess moisture from the chamber. The condensed water tank and the reservoir are both monitored by the system and the touch screen display panel informs the users when the tanks need attention for drainage or refilling.
- **Deflagration Protection** A novel locking mechanism for the door prevents excessive pressure build-up inside the chamber and allows built up pressure to safely escape the chamber, decreasing the chance of deflagration from improper use.



## Airflow and Anti-Condensation Screen

Heated air is guided via an interior glass screen on the door of the NINcha this eliminates any condensation and ensures the user always has a clear view into the chamber. Unlike many industrial climate chambers, the chamber does not use air outlets on the interior walls, and instead creates a homogeneous gentle airstream throughout the entire chamber.




## Integrated Illumination

The development process must be easily observable with a forensic climate chamber. NINcha features both an innovative condensation screen and a large viewing panel in the door. In addition, the chamber has interior illumination on each shelving level that can be turned on as needed



## Touch Panel Display and Process Data Logger

NINcha is microprocessor-controlled by the user via the integrated touch screen display. For ease of use, pre-programmed development cycles for Indandione, Ninhydrin, and DFO can be selected. If required, users can also easily alter the parameters and timer settings via a manual mode for customized cycle options.



The USB data logger records all relevant process data onto a USB pen drive for improved quality control, documentation, and analysis with our included Excel-compatible data analysis software



## Variable Shelving Concept

The brackets on the inner chamber walls of NINcha are built to support both the stainless-steel metal grills and hanging rods for evidence. The grills are made with only two rungs, in order to keep cleaning simple. To support heavier evidence, hanging rods can be used in combination with the grills. An additional bracket at the top of the chamber allows access to the chamber's full height.



## Water Supply

NINcha features an internal tank for the distilled water required for the humidification process. Condensed water produced by the process is collected in a removable container at the bottom of the chamber.



CONTACT OPIRA NOW  
**1300 157 969**  
SOLUTIONS@OPIRA.COM.AU

OPIRA BRISBANE  
32 DIVIDEND ST  
MANSFIELD, QLD 4122

OPIRA MELBOURNE  
25 GRAHAM RD  
CLAYTON SOUTH, VIC 3169

OPIRA PERTH  
9 EARLSTON PLACE  
BOORAGOON, WA 6154

OPIRA SYDNEY  
20 DUKE ST  
FORESTVILLE, NSW 2087



# NINcha

The NINcha series is a family of forensic climate chambers, focusing particularly on the development of fingerprints on porous surfaces after treatment with Ninhydrin, DFO or Indandione.

Evidence like this is usually treated with Ninhydrin, DFO or Indandione by bathing or spraying. Following this chemical treatment process, the samples need to be developed under highly specific temperature and humidity conditions for a set period for optimal fingerprint development results.

NINcha makes all this possible, with a user-friendly climate chamber system, complete with a multilingual touch screen control panel, a USB data logger for quick recording and monitoring of processes, and a novel innovative air stream concept that circulates chemical developer evenly throughout the chamber while avoiding any evidence displacement. Additional features allow users to customize filter settings or use external ventilation systems, customize process settings, and minimize potential damage to evidence. These features both increase user safety and help preserve evidence and prevent evidence contamination.


Other climate chambers on the market operate in a similar fashion to NINcha, but unfortunately often neglect safety aspects and the ability to continuously monitor samples, both of which are extremely important for forensic equipment. NINcha is specially designed for forensic use and built to meet the rigorous requirements of standardized forensic laboratories.

## Technical Data

	NINcha S31	NINcha M31	NINcha L31
<b>External Dimensions (HxWxD)</b>	105 x 70 x 60 cm (with adjustable legs)	185 x 70 x 60 cm (on castors)	205 x 70 x 60 cm (on castors)
<b>Internal Dimensions (HxWxD)</b>	50 x 48 x 48 cm	100 x 48 x 48 cm	150 x 48 x 48 cm
<b>Levels</b>	2 + additional top hanging level	4 + additional top hanging level	6 + additional top hanging level
<b>Temperature Range</b>	25°C <sup>1</sup> - 110° C		
<b>Humidity Range</b>	40 – 80% RH <sup>2</sup> (and humidifier OFF)		
<b>Internal Illumination</b>	2 lamps (one lamp per level)	4 lamps (one lamp per level)	6 lamps (one lamp per level)
<b>Basic Configuration / Optional Features</b>			
<b>Filter System</b> LFD31/LFN31/LFI31	☒/✓/☒	☒/✓/☒	☒/✓/☒
<b>USB Datalogger</b> USB31	✓	✓	✓
<b>UV Decontamination</b> UVC-X1	☒	☒	☒
<b>External Ventilation</b> Upgrade EAA31	✓	✓	✓
<b>Electrical Data</b>			
<b>Voltage:</b>	230V AC/50Hz (110-120V AC or 60Hz version on request)		
<b>Power Requirements:</b>	max. 2.200 W		

**Current:** ca. 10 A (on 230V AC)

- <sup>1</sup> min. temperature is ambient temperature + approx. 5° ✓ fitted in basic configuration  
<sup>2</sup> RH range depending on selected temperature or ○ available option, additional charge



NINcha S31 is a desktop system with adjustable legs (above). An optional storage box BOX31 for any accessories not being used (below, left) is also available.

For each cabinet we offer a UV decontamination unit adapted in power output to the cabinet size (below, right, e.g. UVC-X1 for NINcha M31).

