

## Passbox



### TECHNICAL SPECIFICATIONS

| Model                    | OPIRA-SPB-01                 | OPIRA-SPB-02 | OPIRA-SPB-03 | OPIRA-PB-01         | OPIRA-PB-02 | OPIRA-PB-03 |
|--------------------------|------------------------------|--------------|--------------|---------------------|-------------|-------------|
| External size (W*D*H) mm | 660*570*1000                 | 760*670*1100 | 860*770*1200 | 660*570*580         | 760*670*680 | 860*770*780 |
| Internal size (W*D*H) mm | 500*500*500                  | 600*600*600  | 700*700*700  | 500*500*500         | 600*600*600 | 700*700*700 |
| Material                 | 304 stainless steel          |              |              | 304 stainless steel |             |             |
| Interlock                | Electronic                   |              |              | Electronic          |             |             |
| UV lamp                  | 30W*1                        |              |              | 30W*1               |             |             |
| HEPA filter              | Efficiency≥ 99.999% at 0.3μm |              |              | Not applicable      |             |             |
| Clean air speed          | 0.36~0.54m/s                 |              |              | Not applicable      |             |             |
| Gross weight             | 75kg                         | 100kg        | 130kg        | 35kg                | 40kg        | 50kg        |
| Package size             | 700*770*1200                 | 800*870*1300 | 900*970*1400 | 700*700*780         | 800*800*880 | 900*900*980 |



Our products have been designed and tested to comply with all applicable ACMA regulatory requirements and bear the RCM mark.