

A R A Handy®

Pipette Aid



Compatible with most
0.1-100 mL plastic
and glass pipettes



Li-ion Battery: 2-3 hour
charge for 8+ hours of use.

Lightweight,
ergonomic design
for effortless one-
handed use.

Weight : **200 g**



Autoclaveable
silicone adaptor



**0.45 µm
hydrophobic
filter**

Aspiration

Dispensing

8-step adjustable
aspiration and
dispensing speeds.

Powerful pump:
Aspirates 25 mL
in 5 seconds.

Large LCD Display



- Easily monitor low battery warning and speed settings.
- Usable while charging.

QUALITY | RELIABILITY | VALUE

ARA Handy®

Pipette Aid



ARA Handy Pipette Aid Specifications

Aspiration Speed	25 mL < 5s (8 shift)
Dispensing Speed	Motor (8 shift) / Gravity
Battery	Lithium-ion (3.7 V, 1100mA)
Battery Service Life	More than 8 Hours of Intermittence Use
Charging Time	2-3 Hours
Pipette Types	Glass or Plastic Pipette (0.1-100mL), Pasteur Pipettes
Filter	0.45 µm Hydrophobic
Power Supply / Consumption	Input AC 100-240V / Output voltage 4.2V, 500mA
Weight	200g

Ordering Information

Cat. No.	Description
APD-AID	ARA Handy Pipette Aid (220V, 50/60Hz)
APD-AID-17000103	0.45um Filter
APD-AID-17000205	Lithium battery
APD-AID-17000207	Silicone adapter
APD-AID-17100379	Nose cone

Components



ARA Handy®
Pipette Aid



AC-adaptor
Input AC 100-240V /
Output voltage 4.2V, 500 mA



0.45µm
hydrophobic filter



Wall-mounted Holder
& Two-position Stand



User Manual



Quality Control
Certificate