

Pathology Workstation

Our pathology sampling bench is widely used in hospital pathology departments, laboratories, etc. The ventilation system protects operators from gas produced by the formalin during sampling. In addition hot and cold water system ensures it can adapt to work in different climates.

Our range of pathology workstations features

- Designed and manufactured to comply with AS 2252.6-2011 Controlled environments Clean workstations - Design, install and use
- Adaptable louvre suction inlet - ensures above blower suction and avoids air back-streaming
- Automatic flushing device, slightly right-side sloped worktable for easy waste discharge and cleaning
- Odour-removing function
- Single water sink, two-person operation
- Exhaust air system: High-power suction device, low noise
- An optional water heater is available for hot & cold water adjustment.



Our products have been designed and tested to comply with all applicable ACMA regulatory requirements and bear the RCM mark.

Technical Specifications

MODELS	OPIRA-PWS-1000	OPIRA-PWS-1500	OPIRA-PWS-1800
External Size (W*D*H)	1000 x 750 x 1950 mm	1500 x 750 x 1950 mm	1800 x 750 x 1950 mm
Working Area Size (W*D*H)	895 x 430 x 900 mm	1395 x 430 x 900 mm	1695 x 430 x 900 mm
Work Surface Height	900mm		
Consumption	250W	300W	300W
Optional Water Heater	1500W		
Power Supply	AC220V±10%, 50/60Hz 110V±10%, 60Hz		
Exhaust Airflow Volume	500~600m³/h		
Material	Main Body	304 stainless steel	
	Work Table	304 stainless steel	
Standard Accessories	Spotlight, UV lamp, Fluorescent lamp, exhaust system, 4m exhaust duct, pipe strap, auto flusing device, stainless steel hanger, adjustable hot and cold water tap, 180mm deep stainless steel water sink, funnel with filter, active carbon filter		
Optional Accessories	Pulveriser, electrical water heater, automatic induction hand dryer, eye washer, camera, flex magnifier light, laptop bracket		
Gross Weight	250kg	300kg	340kg
Package size (W*D*H)	1160 x 890 x 2170mm	1660 x 890 x 2170mm	1960 x 890 x 2170mm



INTERNAL CAMERA



FLEXI MAGNIFIER LIGHT (OPTIONAL)



UV LIGHT CONFIGURATION



CONTACT US

1300 157 969

SOLUTIONS@OPIRA.COM.AU

OPIRA BRISBANE

32 DIVIDEND STREET
MANSFIELD, QLD 4122

OPIRA PERTH

9 EARLSTON PLACE
BOORAGOON, WA 6154

OPIRA MELBOURNE

25 GRAHAM ROAD
CLAYTON SOUTH, VIC 3169

OPIRA SYDNEY

20 DUKE ST
FORESTVILLE, NSW 2087