

CSC Chemical Storage Cabinet

A **secure, filtered**, chemical
storage cupboard

FUME [] CARE



CONTACT OPIRA NOW
1300 157 969
SOLUTIONS@OPIRA.COM.AU

HEAD OFFICE - BRISBANE
32 DIVIDEND ST
MANSFIELD, QLD 4122

MELBOURNE
PERTH
SYDNEY



CSC Chemical Storage Cabinet

The CSC range of benchtop or freestanding, filtered storage cupboards are ideal for use in laboratory environments where harmful substances are stored.

Harmful fumes pass through the unit's filters, re-circulating clean air back into the laboratory. There is no need for expensive ducting, as the unit uses a recirculatory airflow. Its ductless construction allows the unit to be easily moved.



Harmful fumes pass through the units filters, re-circulating clean air back into the laboratory.



Carbon and HEPA filters are available to suit most laboratory applications.



There is no need for expensive ducting to be installed due to the units recirculatory airflow.



Its ductless construction allows for the unit to be easily moved.



Specifications

Catalogue No	FUMCSC-86S	FUMCSC-86T
Dimensions (WxDxH), mm	860 x 500 x 820	860 x 500 x 1875
Shelves	1	4
**Main Filter 5kg (qty)	1	1
Airflow (m3/hr)	200	200
Pre-filter (efficiencies superior to 99.6%)	1	1
Electrical Supply	Single Phase, 230v, 50Hz	
Noise (dBA)	55	
Switches	Mains ON/OFF	
Monitoring	Low Airflow Indicator	
Cabinet Construction	Low Noise Centrifugal	
Fan	Epoxy coated steel with lockable doors	
Windows	Clear Acrylic Windows	

Custom Request

Cabinets can also be customised to your specification, please contact us to discuss your requirements.

Consumables

Catalogue No	Description
CSC-PRF	Pack of 6 Prefilters

Filters

Filter Type	GP	FOR	ACI	ETH	SUL	AMM	HEPA	Custom*
CSC-PRF	FCS-001	FCS-013	FCS-015	FCS-018	FCS-016	FCS-012	FCS-030	FCS-110

*Maximum mix of 2 carbons are available on this filter type

To ensure operator protection, the correct filter selection is critical. When ordering please refer to the filter selection chart and also state the contaminants being used.

UK DESIGN AND MANUFACTURE

Designed and built in our Merseyside ISO 9001 accredited facility.

Distributed in Australia by:



1300 157 969
sales@opira.com.au
www.opira.com.au