

Qualification Guide – Opira Laminar Flow Cabinets

General questions

- What is your application?
- What are you currently using?
- When do you need this by/timeline?
- What is your budget?

Specific Questions

- Do your operators require protection from their samples? (if yes, BSC)
- Do your samples require protection/sterility from the environment?
- Are you using any chemicals in the cabinet? (if yes, possibly fume hood)
- What size do you require? (1100, 1500 or 1800 mm)
- Please provide dimensions of the available space (width, depth, ceiling height)
- Do you require a stand? If so, what height? Manual or electric?
- Do you require NATA certification?
- Are there any specific regulations/standards that you must adhere to?
- Would you like us to quote for a service contract/preventative maintenance?

***Site inspection may be required to assess door and ceiling height and available space in lab**

